

Anti-RUNX1T1/ETO/CDR Rabbit pAb



GB114460 -20°C

Product Information			
Protein full name	Protein CBFA2T1		
Synonyms	AML1T1, CBFA2T1, CDR, Cyclin D related protein, Eight twenty one protein, ETO, MTG8, MTG8b, Protein CBFA2T1		
Immunogen	Recombinant protein corresponding to Mouse RUNX1T1		
Uniprot ID	Q61909		
Isotype	IgG		
Purity	Affinity purification		
Predicted MW./Observed MW.	64 kDa / 74 kDa		
Applications	Species	Dilution	Positive Tissue
WB	Mouse, Rat	1: 300-1: 500	brain
Storage			
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Storage Buffer	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.		
background			
Transcriptional corepressor which facilitates transcriptional repression via its association with DNA-binding transcription factors and recruitment of other corepressors and histone-modifying enzymes. Can repress the expression of MMP7 in a ZBTB33-dependent manner. Can repress transactivation mediated by TCF12 (By similarity). Acts as a negative regulator of adipogenesis.			
Images			



WB analysis of RUNX1T1/ETO/CDR(GB114460).
Sample: Protein treated by RIPA Lysis Buffer (G2002).
Blocking buffer: 3% Nonfat dry milk (GC310001) in TBST, RT, 1h.
Primary antibody: 1: 500, 4°C overnight.
Secondary antibody: HRP conjugated Goat Anti-Rabbit IgG (H+L) (GB23303), 1: 5000, RT, 1h.

For Research Use Only!

Ver. No.: V2.0-20260420