

RPMI-1640 Complete Medium

Cat. No.: G4535A

Product Information

Product Name	Cat. No.	Spec.
RPMI-1640 Complete Medium	G4535A-500ML	125 mL×4

Product Description

This product features a formula optimized through long-term testing, which maintains excellent cell growth. It contains all components required for cell culture in the RPMI-1640 system. No additional supplementation is needed; it can be used directly for cell culture. The composition is as follows:

Component	Cat. No.	Volume
RPMI-1640 Medium	G4535	89 % (v/v)
FBS	G8002	10 % (v/v)
Penicillin-Streptomycin Solution (Double Antibody, 100×)	G4003	1 % (v/v)

Storage and Shipping Conditions

Transport in ice packs at low temperature; store at 2–8 °C in the dark. Valid for 6 months.

Note

1. This product is filter-sterilized. Aseptic technique must be practiced during use to avoid contamination.
2. This product contains serum and penicillin-streptomycin solution. No further supplementation of serum or antibiotics is required unless specifically needed.
3. For optimal performance, adhere to the storage conditions. Do not store the product at room temperature or higher for a prolonged period.
4. Use the product within its expiration date. Do not use any product past its expiration date.
5. Do not use if the product appears yellowed or turbid.

For Research Use Only!

Ver. No.: V1.2-20260729